

For Avaira® and Biomedics® XC patients, we recommend upgrading them to:

	Product	Material	Water Content	Base Curve (mm)	Diameter (mm)	Sphere Power	Dk/t (at -3.00D)	Replacement Schedule
1-DAY UPGRADE	MyDay*†	stenfilcon A	54%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	100	1-Day
	clariti® 1 day†	somofilcon A	56%	8.6	14.1	+8.00D to -10.00D (0.50D steps after ±6.00D) No Plano	86	1-Day
2-WEEK UPGRADE	Avaira Vitality™†	fanfilcon A	55%	8.4	14.2	+8.00D to -12.00D (0.50D steps after ±6.00D)	110	2-Week
MONTHLY UPGRADE	Biofinity Energys™	comfilcon A	48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after ±6.00D)	160	Monthly
	Biofinity*†	comfilcon A	48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after ±6.00D)	160	Monthly

For patients in ClearSight™ 1 day toric, Vertex Toric™ XR, Avaira® toric, Frequency® 55 toric and Frequency® 55 toric XR, we recommend upgrading them to:

	Product	Material	Water Content	Base Curve (mm)	Diameter (mm)	Sphere Power	Cylinder Power	Axis	Dk/t (at -3.00D)	Replacement Schedule
1-DAY TORIC UPGRADE	clariti® 1 day toric†	somofilcon A	56%	8.6	14.3	Plano to -9.00D (0.50D steps after -6.00D)	-0.75D -1.25D -1.75D	10° 20° 60° 70° 80° 90° 100° 110° 120° 160° 170° 180°	57	1-Day
						+0.25D to +4.00D (0.25D steps)		20° 70° 90° 110° 160° 180°		
						Plano to -9.00D (0.50D steps after -6.00D)	-2.25D	10° 20° 90° 160° 170° 180°		
2-WEEK TORIC UPGRADE	Avaira Vitality™ toric†	fanfilcon A	55%	8.5	14.5	+8.00D to -10.00D (0.50D steps after ±6.00D)	-0.75D -1.25D -1.75D -2.25D	10° to 180° (10° steps)	90	2-Week
MONTHLY TORIC UPGRADE	Biofinity® toric†	comfilcon A	48%	8.7	14.5	+8.00D to -10.00D (0.50D steps after ±6.00D)	-0.75D -1.25D -1.75D -2.25D	10° to 180° (10° steps)	116	Monthly
	Biofinity® XR toric	comfilcon A	48%	8.7	14.5	+8.50D to +10.00D (0.50D steps)	-0.75D -1.25D -1.75D -2.25D	5° to 180° (5° steps)	116	Monthly
+10.00D to -10.00D (0.50D steps after ±6.00D)	-2.75D -3.25D -3.75D -4.25D -4.75D -5.25D -5.75D									

† And associated private label products.